



DataTherm® II CONTINUOUS TEMPERATURE MONITOR



ACCURACY and FUNCTIONALITY TESTING

Geratherm Medical AG, the manufacturer of the DataTherm® II Continuous Temperature Monitor, specifically recommends that the accuracy of your DataTherm II Monitor be validated every two years.

RG Medical Diagnostics, the Master Distributor of the DataTherm II in the U.S. and Canada, would like to offer our services to test the accuracy and functionality of your DataTherm II monitors.

Our Quality Control staff will conduct the test at our Michigan location following the manufacturer's recommended protocol and accuracy standards using certified testing equipment. Although the DataTherm II was not designed to allow for recalibration, you can feel confident that the tested and approved units you are using in the field are precise and accurate upon completion of this process.

Simply notify us that you will be sending us your DataTherm II monitor for testing...you may call us toll free at 888-596-9498 x203. We will complete the testing within 3-5 business days of receiving the DataTherm and return it to you postage paid. You will receive a signed copy of the test certificate for your records and we will retain a test certificate as well listing your DataTherm by serial number in our files. There is a modest fee of \$65.00 for our service per unit. **(This will include a certification document attesting to the accuracy and functionality of the device).**

Please contact me for any questions you may have and or to schedule accuracy testing for your DataTherm monitor.

Mary LeTourneau, RN

RG Medical Diagnostics

28351 Beck Road, Suite G5

Wixom, Michigan 48393

248-504-6726 ext. 207

maryRN@rgmd.com

