

What is Tecartherapy?

T.E.C.A.R. means Transfer of Capacitive And Resistive Energy. It is a therapy of high frequency current/radiofrequency acting at cells level with **3 key therapeutic effects**:

› Biological


boosts cellular metabolism/permeability

› Analgesic

Pain relief & nociceptors saturation

› Diathermy

superficial/deep heat, long lasting warmth feeling



DEEP ACTION

HARD TISSUE = RET

LIGAMENT, JOINT, TENDON, PERIOSTEUM, BONES

1 MHz

Superficial Action



500 kHz

Semi Deep Action




300 kHz

Deep Action

SUPERFICIAL ACTION

SOFT TISSUE = CET

MUSCLE, FASCIA, APONEUROSIS, VISCERAL, LYMPHATIC AND VASCULAR.



WINBACK TECAR IS A PHYSIOLOGICAL BIO ACCELERATOR

◀ **ELIMINATE PAIN**

Inhibits pain immediately for a duration of more than 48 hours after treatment

◀ **FREE UP MOVEMENT**

Relieves tension by revascularising them locally. Your range of movement is immediately restored

◀ **ACCELERATE HEALING**

Accelerates the body's natural self-healing system by promoting intra and extracellular exchanges in a lasting way

- ◀ It is **non-invasive** and 100% natural for the body
- ◀ It offers a **gentle patient experience** by eliminating pain immediately
- ◀ It **is fully adapted to therapists' techniques**: manual, mobilisation, active
- ◀ Limited contraindications
- ◀ Can be combined with other popular tools or techniques such as **IASTM, Cupping, Taping, ART** also makes other 'passive' modalities better if layered in



TECAR Packs	BACK ¹	BACK ³
TECAR 1.0	✓	✓
TECAR 2.0	✓	✓
TECAR 3.0	✓	✓
TECAR 4.0	✗	✓
TECAR 5.0	✓	✓
Trolley	✗	✓
Pelican Case Available	✓	✓
Unit frequencies / Function / Power		
CET (500 KHz)	✓	✓
Deep CET (300 KHz)	✓	✓
RET (500 KHz)	✓	✓
RET+ (300 KHz)	✗	✓
MIX (1000 KHz) & MIX+ (500 KHz)	✗	✓
Dynamic	✓	✓
Low Pulse & Super Pulse	✓	✓
TIC & Scan	✓	✓
SWAP & BEAT	✗	✓
Size	12.9" x 12.9" 5.9"	14" X 12" 6
Weight	9 lbs.	9 lbs.
Trolley Size	✗	20" x 13" x 39"
Trolley Weight	✗	37 lbs.
Power CET	400 VA	400 VA
Power RET	100 Wats	100 Watts